Proximal Cross Locking Tibial Intramedullary Nails PCL

V2-EVREN

Sanayi Mah. Sancaklı Cad. Başaklı Sok. No:20
34165 Güngören - İstanbul / TÜRKİYE

Tel: +90 212 644 53 42 Pbx.

Fax: +90 212 504 02 28

e-mail: sales@v2evren.com.tr

http://www.v2evren.com.tr

Proximal Cross Locking Tibial Nails (PCL) St. Steel

Length mm	Nail Diameters					
	Ø8mm	Ø9mm	Ø10mm	Ø11mm	Ø12mm	Ø13mm
200 mm	TB-2822					
240 mm	TB-2824	TB-2924	TB-2024	TB-2124		
260 mm	TB-2826	TB-2926	TB-2026	TB-2126	TB-2226	TB-2326
280 mm	TB-2828	TB-2928	TB-2028	TB-2128	TB-2228	TB-2328
300 mm	TB-2830	TB-2930	TB-2030	TB-2130	TB-2230	TB-2330
320 mm	TB-2832	TB-2932	TB-2032	TB-2132	TB-2232	TB-2332
340 mm	TB-2834	TB-2934	TB-2034	TB-2134	TB-2234	TB-2334
360 mm	TB-2836	TB-2936	TB-2036	TB-2136	TB-2236	TB-2336
380 mm	TB-2838	TB-2938	TB-2038	TB-2138	TB-2238	TB-2338
400 mm	TB-2840	TB-2940	TB-2040	TB-2140	TB-2240	TB-2340

Nails in the shaded areas are available on special request.

Proximal Cross Locking Tibial Nails (PCL) Pure Titanium

Length mm	Nail Diameters					
	Ø8mm	Ø9mm	Ø10mm	Ø11mm	Ø12mm	Ø13mm
200 mm	TB-2822T					
240 mm	TB-2824T	TB-2924T	TB-2024T	TB-2124T		
260 mm	TB-2826T	TB-2926T	TB-2026T	TB-2126T	TB-2226T	TB-2326T
280 mm	TB-2828T	TB-2928T	TB-2028T	TB-2128T	TB-2228T	TB-2328T
300 mm	TB-2830T	TB-2930T	TB-2030T	TB-2130T	TB-2230T	TB-2330T
320 mm	TB-2832T	TB-2932T	TB-2032T	TB-2132T	TB-2232T	TB-2332T
340 mm	TB-2834T	TB-2934T	TB-2034T	TB-2134T	TB-2234T	TB-2334T
360 mm	TB-2836T	TB-2936T	TB-2036T	TB-2136T	TB-2236T	TB-2336T
380 mm	TB-2838T	TB-2938T	TB-2038T	TB-2138T	TB-2238T	TB-2338T
400 mm	TB-2840T	TB-2940T	TB-2040T	TB-2140T	TB-2240T	TB-2340T

Nails in the shaded areas are available on special request.

PCL Locking Screws Pure Titanium and Stainless Steel

Cat. No	Cat. No	Lenght	Tread Lenght
TB-1420T	TB-1420	20	7mm
TB-1425T	TB-1425	25	
TB-1430T	TB-1430	30	
TB-1435T	TB-1435	35	
TB-1440T	TB-1440	40	
TB-1445T	TB-1445	45	

Cat. No	Cat. No	Lenght	Tread Lenght
TB-1450T	TB-1420	20	12mm
TB-1455T	TB-1425	25	
TB-1455T	TB-1430	30	
TB-1465T	TB-1435	35	
TB-1470T	TB-1440	40	
TB-1475T	TB-1445	45	

PCL End Cup Pure Titanium and Stainless Steel

TB-1401

TB-1401T

ÇAKMA TAKIMI / INSTRUMENTS

TB-8010 Handle + Bar Locking Screw

TB-8015 Guide Bar

TB-8020 Distal Outrigger

TB-8025 Proximal Outrigger

TB-8035 Locking Rod

SN-8090 Supine Driver

TB-8040 Allen Wrench 5 mm

TB-8120 Sliding Hammer With Detachable Swing arm

TB-8125 Blue Handle With Bayonet Fitting

TB-8075 Allen Wrench 3 mm

TB-8045 Screw Guides 10/8

TB-8050 Drill Guides Ø4 mm

TB-8055 Graduated Angled Trocar

FNI-3430 Drill Bit 4.0 mm

TB-8090 T-Handle Reamer, Ø4 mm

TB-8065 Drill Stopper Ø4 mm

U-Shaped Stabilizing Spacers

TB-8080-1 8/13

TB-8080-2 9/12

TB-8080-3 10/11

TB-8113 Open Wrench, 13 mm

TB-8130 Metal Ruler

TB-8070 Locking Screw Depth Gauge

TB-8085 T-Handled Stabilizing Rod

TB-8095 T-Handled Locking Screw Extractor

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
Proximal Cross Locking Tibial Nail has been designed for Tibial shaft fractures within 4 cm distal to the Tibial tuberosity to 4cm proximal to the ankle joint tuberosity.

Tibia kırıklarında sınıflama (Chapman) Classification of tibial fractures (Chapman)

- A- Transvers / kısa oblik - **Transverse / short oblique**
- B- Küçük Kelebek - **Small butterfly**
- C- Geniş Kelebek - **Large butterfly**
- D- Segmenter - **Segmental**
- E- Spiral - **Spiral**
- G- Distal - **Distal**

Tibia Diafiz kırıklarında parçalanma sınıflaması
(Henley Mb Intramedullary devices for Tibial fracture stabilisation Clin Orthop. 240:87-97-1989) **Classification of Tibial shaft comminution**

Kırıkta kortikal temas en az %50 ve 1/3 diafiz orta hatta ise stabil olarak kabul edilir.
Fracture is stable with at least 50% of cortical continuity and within the middle 1/3 of the Tibia.

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When stability is in doubt static interlocking should be the choice to control the alignment, rotation and length.

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Measuring length of nail is performed with a tape measuring length from medial malleolus to the Tibial tuberosity on the unaffected site nail should come within 2cm of articular surface at the end of tibia to provide maximum fixation

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To check the rotational alignment one should compare with the unaffected leg.

Cerrahi teknik Surgical Technique

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
Traction is applied either with a calcaneal traction pin or holding the foot.

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Be sure not to put too much pressure on the Popliteal region.

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
A 5 cm incision is made medial to the patellar ligament.

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The patellar tendon can be retracted medially or it can be splitted in the midline.

Pain and heterotopic ossification are reported after midline splitting technique.

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In older patients or with proximal fractures the original entrance site can cause posterior blowout. A more proximal portal beginning approximately 3-4 mm back on the superior surface of tibia is closer to the intramedullary line and avoids these problems.

PCL TİBİA ÇİVİ SİSTEMİ AMELİYAT TEKNİĞİ / PCL TIBIAL NAILING SYSTEM SURGICAL PROTOCOL

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There are two entry points: Superior and anterior. The preferred entry point is the superior approach since it allows easier alignment with the medullary canal. A hole is made with TB-8115 awl. When this approach is used the entry portal must be very proximal no more than 1cm distal to the anterior edge of the tibial plateau. A more distal entry point may result in damage to the posterior cortex. A hole is made with TB-8115 awl.

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Push the TB-8135 guide wire until its tip sits 0.5-1cm proximal to the ankle joint (talotibial) taking care to ensure that it is exactly in the midline.

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Ream to a width 1-2 mm greater than the proposed nail (use the FNI-3007, 3008, 3009, and 3010 flexible reamers.)

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Insert the TB-8035 Locking Rod into the TB-8010 back of the handle and the nail of correct diameter and length into the nail support and tighten the locking rod with the TB-8040 5mm Allen wrench.

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Insert the nail, if reamed over the guide wire, as far as possible using image intensification. Do not forget that there is no guide hole in 8 and 9 mm diameter nails. Remove guide wire, if applicable, when its exit point from the nail is at the level of the entry portal. The nail is correctly inserted when the step of the nail support is flush with the surface of the bone.

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The TB-8120 Sliding Hammer, attached to the end of the nail locking rod end fully tightened, may be used to insert the nail gently in the correct position. Check that the TB-8035 locking rod is tight after hammer is removed.

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Insert the TB-8015 guide bar into the handle, adjust its position until the number corresponding to the selected nail length lines up with the front of the handle. Lock the guide bar firmly into place.

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Mount the TB-8020 distal outrigger on the guide bar so that it lies on the correct side of the tibia, usually medial. Insert the TB-8045 10/8 screw guides into the outrigger. No incision is made yet!

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Insert the TB-8050 drill guide into the hole in the guide bar just proximal to the distal outrigger. Make an incision and advance the drill guide until its teeth are engaged in the tibia and stabilized on the center of the tibial crest. A 4mm drill bit (TB-8060) is used to drill the anterior cortex only.

PCL TIBIA ÇİVİ SİSTEMİ AMELİYAT TEKNİĞİ / PCL TIBIAL NAILING SYSTEM SURGICAL PROTOCOL

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
Clear the hole in the bone with the TB-8090 4 mm T-handled Reamer until the reamer can be heard tapping the nail.

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Remove the TB-8090 4 mm T-handled Reamer and insert the TB-8085 T-handled Stabilizing Rod, down to the nail, again tapping the nail to ensure contact.

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Maintain contact between the tip of the Stabilizing Rod and the nail. The stabilizing rod may have to be lifted up or pushed down to establish correct contact with the nail. Make an incision beneath each screw guide. Advance the screw guides until they are in contact with the cortex. Tighten the clamp locking nut on the outrigger to hold them firmly in place. Insert the 4mm drill guide into the most distal of the screw guides and drill the bone with the 4mm drill bit. While the surgeon is drilling, the assistant must hold the T-handle of the Stabilizing Rod, keep its tip against the nail, and maintain this position throughout the drilling procedure.

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Remove the TB-8050 drill bit and drill guide, and insert the TB-8055 8/4 mm graduated angled trocar. Drill the second hole in the same way.

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Insert the locking screw of correct length. If you need to take out the screws you may use the TB-8095 T Handled screw remover. Be sure to take out the screw guide first. Use the TB-8100 T Handled 3.5mm screw- driver. Remove the distal outrigger and the T-handled stabilizing rod.

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Check for any malrotation or distraction of the fracture site, before carrying out proximal locking. If the fracture site is distracted, attach the TB- 8120 sliding hammer to the locking rod and close the fracture gap by gentle reverse hammering.

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Attach the correct Stabilizing Spacer (TB-8080-1,-2,-3) for the diameter of the nail to the T-handled Stabilizing Rod. Position the Spacer so that the correct nail diameter is visible on the upper surface, facing towards the surgeon.

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Loosen the guide bar locking screw and move the guide bar until the "P" mark is level with the front surface of the handle. Lock the bar into position.

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Mount the TB-8025 proximal outrigger on the guide bar and insert two screw guides into the holes in the proximal outrigger. Make an incision and advance the screw guides down to the cortex. Lock them into position with the clamp locking nuts.

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The medial hole is drilled first. Insert the 4 mm drill guide into the screw guide and drill the bone with the 4 mm drill bit. After drilling, remove the drill bit and drill guide, and insert the TB-8055 graduated angled trocar. Drill the lateral hole and insert the locking screws of correct length.

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Remove the proximal outrigger and the guide bar. Before removing the handle from the nail, check correct insertion of locking screws both in AP and lateral planes. Remove the TB-8035 locking rod and the handle (use the TB-8040 AA 5 mm Allen wrench if necessary.)

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Take out the TB-1401 nail end cap first with the TB-8100 screwdriver. Mount the adaptor manually to the nail. If there is a need for more stability than you may use the TB-8113 wrench.

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
The locking screws can be taken out proximally and distally with TB-8095 and TB-8100 screw drivers.

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The extractor handle TB-8120 is attached to the screw adaptor, and the nail is then removed from the tibia.

V2-EVREN

Tibbi Mamüller İmalat ve Ticaret A.Ş.

Sanayi Mah. Sancaklı Cd. Başaklı Sk.
No 20 34610 Güngören İstanbul / TÜRKİYE
Tel:+90 212 644 53 42 Pbx.
Fax:+90 212 504 02 28
e-mail:sales@v2evren.com.tr
<http://www.v2evren.com.tr>